



**BEESWIFT®**  
FOCUSED ON SAFETY

Product Code  
CM1650

Available in Sizes  
**ONE SIZE**

Available in Colours  
 White  
CM1650

# DATASHEET

## Zoll AED 3 CPR Adult/Child Uni-Padz




- For use exclusively with the ZOLL AED 3®, CPR Uni-padz electrodes have an industry-leading shelf life of five years and can be used on both adult and child patients, so there is no need to maintain two separate sets of pads.
- The ZOLL AED 3 full-colour LCD display includes graphics that coincide with the electrode shapes and colours to support proper pad placement on both adults and children.
- By simply repositioning the pads and activating the Child Mode on the ZOLL AED 3, CPR Uni-padz can be used to treat child patients 8 years of age or younger or under 25 kg (55 lbs.).

### CPR Uni-padz™


1



2



3



OPERATING TEMPERATURE: 0°C to 50°C (32°F to 122°F)  
SHORT TERM STORAGE TEMPERATURE: -30°C to 65°C (-22°F to 149°F)  
LONG TERM STORAGE TEMPERATURE: 0°C to 35°C (32°F to 95°F)

 Rx ONLY

 0°C to 50°C  
32°F to 122°F

 NON-LATEX

 NO FIRE

 NO SHOCK

 NO WATER

 CE 0123

R2138-30 Rev. B Patent: [www.zoll.com/patents](http://www.zoll.com/patents)



ZOLL Medical Corporation  
2400 AED Road  
Chatham, MA 01924-1001 USA  
800.348.9011  
Made in USA

ZOLL International Holding B.V.  
Neukirchweg 18 6662 PV EIJH  
The Netherlands Tel: +31 (0) 481 366410

LOT 2321

REF 8900-000270

2026-09-05

(01)00847946024161(17)260905(10)2321



